

# Certificate of performance



**Article:** 3 715 001

**Grad. cylinder, 10:0,2 ml, A, conf.cert.**




**Batch:** 06-18

**Serial:**

**Measuring Nr.** 17998

**Batch size:** 1000

Volume:	2,0000 ml	6,0000 ml	10,0000 ml
Tolerance:	± 0,1000 ml	± 0,1000 ml	± 0,1000 ml
Average value:	1,9925 ml	6,0180 ml	10,0331 ml
Deviation:	0,0165 ml	0,0096 ml	0,0137 ml
Precision (CV):	0,8272 %	0,1594 %	0,1370 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

Operator:	9
Date:	14.02.2018
Balance:	Sartorius - LE244S-OCE

