

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 715 001

**Grad. cylinder, 10:0,2 ml, A, conf.cert.**




**Batch:** 10-20

**Serial:**

**Measuring Nr.** 36489

**Batch size:** 1200

Volume:	2,0000 ml	6,0000 ml	10,0000 ml
Tolerance:	± 0,1000 ml	± 0,1000 ml	± 0,1000 ml
Average value:	1,9949 ml	6,0167 ml	10,0344 ml
Deviation:	0,0203 ml	0,0193 ml	0,0213 ml
Precision (CV):	1,0167 %	0,3216 %	0,2135 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

Operator:	9
Date:	05.03.2020
Balance:	Sartorius - LE244S-OCE

