

Certificate of performance



Article: 3 715 001

Grad. cylinder, 10:0,2 ml, A, conf.cert.




Batch: 27-21

Serial:

Measuring Nr. 47786

Batch size: 1200

Volume:	2,0000 ml	6,0000 ml	10,0000 ml
Tolerance:	± 0,1000 ml	± 0,1000 ml	± 0,1000 ml
Average value:	1,9992 ml	6,0021 ml	10,0013 ml
Deviation:	0,0172 ml	0,0157 ml	0,0182 ml
Precision (CV):	0,8579 %	0,2612 %	0,1823 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	08.07.2021
Balance:	Sartorius - LE244S-OCE

