

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 003 BL USP

Grad. cyl., 25:0,5 ml, blue grad., USP, A, conf. cert.




Batch: 22-20

Serial:

Measuring Nr. 38353

Batch size: 150

Volume:	5,0000 ml	15,0000 ml	25,0000 ml
Tolerance:	± 0,1700 ml	± 0,1700 ml	± 0,1700 ml
Average value:	4,9810 ml	15,0291 ml	25,0543 ml
Deviation:	0,0236 ml	0,0344 ml	0,0367 ml
Precision (CV):	0,4720 %	0,2295 %	0,1467 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	02.06.2020
Balance:	Sartorius - MSE224S-0CE -DU

