

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 003 USP

Grad. cyl., 25 : 0,5 ml, USP, A, conf. cert




Batch: 34-23

Serial:

Measuring Nr. 64155

Batch size: 150

Volume:	5,0000 ml	15,0000 ml	25,0000 ml
Tolerance:	± 0,1700 ml	± 0,1700 ml	± 0,1700 ml
Average value:	5,0054 ml	14,9926 ml	24,9822 ml
Deviation:	0,0289 ml	0,0158 ml	0,0144 ml
Precision (CV):	0,5777 %	0,1051 %	0,0576 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	26.08.2023
Balance:	Sartorius - MSE224S-0CE -DU

