

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 715 003

**Graduated Cylinder, 25:0,5 ml, A, conf. cert**




**Batch:** 28-17

**Serial:**

**Measuring Nr.** 12118

**Batch size:** 800

| Volume:         | 5,0000 ml   | 15,0000 ml  | 25,0000 ml  |
|-----------------|-------------|-------------|-------------|
| Tolerance:      | ± 0,2500 ml | ± 0,2500 ml | ± 0,2500 ml |
| Average value:  | 5,0432 ml   | 15,0602 ml  | 25,0639 ml  |
| Deviation:      | 0,0715 ml   | 0,0632 ml   | 0,0576 ml   |
| Precision (CV): | 1,4297 %    | 0,4211 %    | 0,2303 %    |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

|           |                           |
|-----------|---------------------------|
| Operator: | 44                        |
| Date:     | 19.07.2017                |
| Balance:  | Sartorius -<br>LE244S-OCE |

