

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 715 003

**Graduated cylinder, 25:0,5 ml, A, conf. cert**




**Batch:** 22-20

**Serial:**

**Measuring Nr.** 38463

**Batch size:** 1200

Volume:	5,0000 ml	15,0000 ml	25,0000 ml
Tolerance:	± 0,2500 ml	± 0,2500 ml	± 0,2500 ml
Average value:	5,0131 ml	15,0306 ml	25,0385 ml
Deviation:	0,0506 ml	0,0452 ml	0,0484 ml
Precision (CV):	1,0112 %	0,3012 %	0,1937 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

Operator:	9
Date:	08.06.2020
Balance:	Sartorius - LE244S-OCE

