

# Certificate of performance



**Article:** 3 715 003

**Graduated cylinder, 25:0,5 ml, A, conf. cert**




**Batch:** 31-22

**Serial:**

**Measuring Nr.** 57819

**Batch size:** 279

Volume:	5,0000 ml	15,0000 ml	25,0000 ml
Tolerance:	± 0,2500 ml	± 0,2500 ml	± 0,2500 ml
Average value:	5,0234 ml	15,0464 ml	25,0738 ml
Deviation:	0,0237 ml	0,0159 ml	0,0166 ml
Precision (CV):	0,4741 %	0,1061 %	0,0663 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

Operator:	7
Date:	04.08.2022
Balance:	Sartorius - LE244S-OCE

