

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 005 BL USP

Grad. cylinder, 50:1 ml, blue grad., USP, A, conf. cert.




Batch: 09-19

Serial:

Measuring Nr.: 28173

Batch size: 50

Volume:	10,000 ml	30,000 ml	50,000 ml
Tolerance:	± 0,2500 ml	± 0,2500 ml	± 0,2500 ml
Average value:	9,9612 ml	30,0000 ml	50,0215 ml
Deviation:	0,0413 ml	0,0374 ml	0,0607 ml
Precision (CV):	0,4127 %	0,1246 %	0,1214 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	9
Date:	07.03.2019
Balance:	Sartorius - MSE224S-0CE -DU

