

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 715 005 BL

**Graduated cylinders, 50:1 ml, blue grad., A, conf. cert.**




**Batch:** 29-18

**Serial:**

**Measuring Nr.** 22424

**Batch size:** 500

| Volume:                | 10,000 ml   | 30,000 ml   | 50,000 ml   |
|------------------------|-------------|-------------|-------------|
| <b>Tolerance:</b>      | ± 0,5000 ml | ± 0,5000 ml | ± 0,5000 ml |
| <b>Average value:</b>  | 9,8801 ml   | 29,9666 ml  | 49,9968 ml  |
| <b>Deviation:</b>      | 0,0856 ml   | 0,0870 ml   | 0,0971 ml   |
| <b>Precision (CV):</b> | 0,8555 %    | 0,2900 %    | 0,1941 %    |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

**Operator:** 44

**Date:** 25.07.2018

**Balance:** Sartorius -  
LE244S-OCE

DE-M