

Certificate of performance



Article: 3 715 005 BL

Graduated cylinders, 50:1 ml, blue grad., A, conf. cert.




Batch: 41-19

Serial:

Measuring Nr.: 33639

Batch size: 280

Volume:	10,000 ml	30,000 ml	50,000 ml
Tolerance:	± 0,5000 ml	± 0,5000 ml	± 0,5000 ml
Average value:	9,9606 ml	30,0045 ml	50,0026 ml
Deviation:	0,0783 ml	0,0742 ml	0,0608 ml
Precision (CV):	0,7833 %	0,2474 %	0,1216 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	14.10.2019
Balance:	Sartorius - LE244S-OCE

