

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 005 BL

Graduated cylinders, 50:1 ml, blue grad., A, conf. cert.




Batch: 40-20

Serial:

Measuring Nr.: 40642

Batch size: 280

Volume:	10,000 ml	30,000 ml	50,000 ml
Tolerance:	± 0,5000 ml	± 0,5000 ml	± 0,5000 ml
Average value:	9,9633 ml	29,9863 ml	49,9974 ml
Deviation:	0,0592 ml	0,0458 ml	0,0443 ml
Precision (CV):	0,5922 %	0,1527 %	0,0885 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	29.09.2020
Balance:	Sartorius - LE244S-OCE

