

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 715 005 USP

**Grad.cylinder, 50:1 ml, USP, A, conf.cert.**




**Batch:** 23-22

**Serial:**

**Measuring Nr.:** 56536

**Batch size:** 280

Volume:	10,000 ml	30,000 ml	50,000 ml
Tolerance:	± 0,2500 ml	± 0,2500 ml	± 0,2500 ml
Average value:	10,0190 ml	30,0260 ml	50,0012 ml
Deviation:	0,0311 ml	0,0347 ml	0,0400 ml
Precision (CV):	0,3105 %	0,1158 %	0,0801 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

Operator:	7
Date:	10.06.2022
Balance:	Sartorius - LE244S-OCE

