

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 005 USP

Graduated cylinder, 50:1 ml, USP, A, conf.cert.




Batch: 23-23

Serial:

Measuring Nr.: 63300

Batch size: 300

Volume:	10,0000 ml	30,0000 ml	50,0000 ml
Tolerance:	± 0,2500 ml	± 0,2500 ml	± 0,2500 ml
Average value:	10,0222 ml	30,0219 ml	49,9667 ml
Deviation:	0,0406 ml	0,0425 ml	0,0416 ml
Precision (CV):	0,4062 %	0,1418 %	0,0832 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	13.06.2023
Balance:	Sartorius - LE244S-OCE

