

Certificate of performance



Article: 3 715 005

Graduated cylinders, A, conf. cert.




Batch: 44-17

Serial:

Measuring Nr.: 15696

Batch size: 500

Volume:	10,000 ml	30,000 ml	50,000 ml
Tolerance:	± 0,5000 ml	± 0,5000 ml	± 0,5000 ml
Average value:	9,9997 ml	30,0678 ml	50,0937 ml
Deviation:	0,0640 ml	0,0381 ml	0,0233 ml
Precision (CV):	0,6395 %	0,1270 %	0,0466 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	13.11.2017
Balance:	Sartorius - LE244S-OCE

