

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 715 005

**Graduated cylinders, 50:1 ml, A, conf. cert.**




**Batch:** 08-20

**Serial:**

**Measuring Nr.** 36403

**Batch size:** 1200

Volume:	10,000 ml	30,000 ml	50,000 ml
Tolerance:	± 0,5000 ml	± 0,5000 ml	± 0,5000 ml
Average value:	10,0579 ml	30,0690 ml	50,0755 ml
Deviation:	0,0597 ml	0,0588 ml	0,0636 ml
Precision (CV):	0,5966 %	0,1962 %	0,1273 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

Operator:	7
Date:	26.02.2020
Balance:	Sartorius - MSE224S-0CE -DU

