

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 005

Graduated cylinders, 50:1 ml, A, conf. cert.




Batch: 40-20

Serial:

Measuring Nr. 40675

Batch size: 500

Volume:	10,0000 ml	30,0000 ml	50,0000 ml
Tolerance:	± 0,5000 ml	± 0,5000 ml	± 0,5000 ml
Average value:	9,9813 ml	29,9869 ml	49,9837 ml
Deviation:	0,0769 ml	0,0506 ml	0,0509 ml
Precision (CV):	0,7694 %	0,1685 %	0,1019 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	30.09.2020
Balance:	Sartorius - LE244S-OCE

