

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 005

Graduated cylinders, 50:1 ml, A, conf. cert.




Batch: 51-21

Serial:

Measuring Nr. 51290

Batch size: 499

Volume:	10,000 ml	30,000 ml	50,000 ml
Tolerance:	± 0,5000 ml	± 0,5000 ml	± 0,5000 ml
Average value:	10,0115 ml	30,0348 ml	50,0522 ml
Deviation:	0,0533 ml	0,0501 ml	0,0363 ml
Precision (CV):	0,5330 %	0,1669 %	0,0726 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	04.01.2022
Balance:	Sartorius - LE244S-OCE

