

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 715 010 BL USP

**Grad. cylinder, 100:1 ml, blue grad. USP, A, conf. cert.**




**Batch:** 11-19

**Serial:**

**Measuring Nr.:** 28317

**Batch size:** 90

Volume:	20,000 ml	60,000 ml	100,000 ml
<b>Tolerance:</b>	± 0,5000 ml	± 0,5000 ml	± 0,5000 ml
<b>Average value:</b>	19,8186 ml	59,8360 ml	99,8393 ml
<b>Deviation:</b>	0,0952 ml	0,0642 ml	0,1117 ml
<b>Precision (CV):</b>	0,4762 %	0,1069 %	0,1117 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

<b>Operator:</b>	0
<b>Date:</b>	15.03.2019
<b>Balance:</b>	Sartorius - LE244S-OCE

