

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 715 010 BL

**Grad. cylinder, 100:1,0 ml, blue grad., A, conf. cert.**




**Batch:** 25-18

**Serial:**

**Measuring Nr.:** 21855

**Batch size:** 700

Volume:	20,000 ml	60,000 ml	100,000 ml
<b>Tolerance:</b>	± 0,5000 ml	± 0,5000 ml	± 0,5000 ml
<b>Average value:</b>	19,9335 ml	59,9654 ml	99,9693 ml
<b>Deviation:</b>	0,0891 ml	0,0512 ml	0,0546 ml
<b>Precision (CV):</b>	0,4453 %	0,0854 %	0,0546 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

**Operator:** 44

**Date:** 02.07.2018

**Balance:** Sartorius -  
LE244S-OCE

