

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 010 BL

Grad. cylinder, 100:1,0 ml, blue grad., A, conf. cert.




Batch: 29-20

Serial:

Measuring Nr. 39195

Batch size: 280

Volume:	20,0000 ml	60,0000 ml	100,0000 ml
Tolerance:	± 0,5000 ml	± 0,5000 ml	± 0,5000 ml
Average value:	19,8968 ml	60,0163 ml	100,1083 ml
Deviation:	0,0847 ml	0,0712 ml	0,0759 ml
Precision (CV):	0,4233 %	0,1186 %	0,0759 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	24.07.2020
Balance:	Sartorius - MSE224S-0CE -DU

