

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 715 010 BL

**Grad. cylinder, 100:1,0 ml, blue grad., A, conf. cert.**




**Batch:** 13-21

**Serial:**

**Measuring Nr.:** 45477

**Batch size:** 500

Volume:	20,000 ml	60,000 ml	100,000 ml
Tolerance:	± 0,5000 ml	± 0,5000 ml	± 0,5000 ml
Average value:	20,0253 ml	60,0233 ml	100,0536 ml
Deviation:	0,0866 ml	0,1074 ml	0,0755 ml
Precision (CV):	0,4329 %	0,1790 %	0,0755 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

Operator:	7
Date:	06.04.2021
Balance:	Sartorius - LE244S-OCE

