

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 010 BL

Grad. cylinder, 100:1,0 ml, blue grad., A, conf. cert.




Batch: 29-23

Serial:

Measuring Nr. 63829

Batch size: 620

Volume:	20,0000 ml	60,0000 ml	100,0000 ml
Tolerance:	± 0,5000 ml	± 0,5000 ml	± 0,5000 ml
Average value:	19,9469 ml	60,0198 ml	100,0002 ml
Deviation:	0,0935 ml	0,1095 ml	0,0969 ml
Precision (CV):	0,4673 %	0,1825 %	0,0969 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	20.07.2023
Balance:	Sartorius - LE244S-OCE

