

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 010 BL

Grad. cylinder, 100:1,0 ml, blue grad., A, conf. cert.




Batch: 47-23

Serial:

Measuring Nr.: 65414

Batch size: 500

| Volume: | 20,0000 ml | 60,0000 ml | 100,0000 ml |
|------------------------|-------------|-------------|-------------|
| Tolerance: | ± 0,5000 ml | ± 0,5000 ml | ± 0,5000 ml |
| Average value: | 19,9334 ml | 59,9777 ml | 100,0608 ml |
| Deviation: | 0,0793 ml | 0,0792 ml | 0,0789 ml |
| Precision (CV): | 0,3963 % | 0,1320 % | 0,0789 % |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

| | |
|------------------|---------------------------|
| Operator: | 7 |
| Date: | 01.12.2023 |
| Balance: | Sartorius - LE244S-OCE |

