

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 715 010 USP

**Graduated cylinders, 100:1 ml, USP, A, conf. cert.**




**Batch:** 25-18

**Serial:**

**Measuring Nr.** 21873

**Batch size:** 70

Volume:	20,000 ml	60,000 ml	100,000 ml
<b>Tolerance:</b>	± 0,5000 ml	± 0,5000 ml	± 0,5000 ml
<b>Average value:</b>	19,9390 ml	59,9349 ml	99,9195 ml
<b>Deviation:</b>	0,1127 ml	0,0739 ml	0,0529 ml
<b>Precision (CV):</b>	0,5637 %	0,1232 %	0,0529 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

<b>Operator:</b>	44
<b>Date:</b>	03.07.2018
<b>Balance:</b>	Sartorius - LE244S-OCE

