

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 010 USP

Graduated cylinders, 100:1 ml, USP, A, conf. cert.




Batch: 09-23

Serial:

Measuring Nr.: 61881

Batch size: 500

Volume:	20,000 ml	60,000 ml	100,000 ml
Tolerance:	± 0,5000 ml	± 0,5000 ml	± 0,5000 ml
Average value:	20,0013 ml	59,9601 ml	99,9406 ml
Deviation:	0,0955 ml	0,0777 ml	0,0645 ml
Precision (CV):	0,4774 %	0,1296 %	0,0645 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	07.03.2023
Balance:	Sartorius - LA 1200 S

DE-M