

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 715 010

**Graduated cylinder, 100 ml, A, conf. cert.**




**Batch:** 29-17

**Serial:**

**Measuring Nr.** 12188

**Batch size:** 1200

Volume:	20,000 ml	60,000 ml	100,000 ml
<b>Tolerance:</b>	± 0,5000 ml	± 0,5000 ml	± 0,5000 ml
<b>Average value:</b>	19,9904 ml	60,0582 ml	100,1151 ml
<b>Deviation:</b>	0,1321 ml	0,1079 ml	0,1073 ml
<b>Precision (CV):</b>	0,6606 %	0,1798 %	0,1073 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

<b>Operator:</b>	9
<b>Date:</b>	24.07.2017
<b>Balance:</b>	Sartorius - LE244S-OCE

