

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 715 010

**Graduated cylinder, 100:1,0 ml, A, conf. cert.**




**Batch:** 10-18

**Serial:**

**Measuring Nr.** 19063

**Batch size:** 500

| Volume:                | 20,000 ml   | 60,000 ml   | 100,000 ml  |
|------------------------|-------------|-------------|-------------|
| <b>Tolerance:</b>      | ± 0,5000 ml | ± 0,5000 ml | ± 0,5000 ml |
| <b>Average value:</b>  | 19,9573 ml  | 59,9677 ml  | 99,9670 ml  |
| <b>Deviation:</b>      | 0,0978 ml   | 0,0820 ml   | 0,0811 ml   |
| <b>Precision (CV):</b> | 0,4889 %    | 0,1367 %    | 0,0811 %    |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

|                  |                           |
|------------------|---------------------------|
| <b>Operator:</b> | 44                        |
| <b>Date:</b>     | 23.03.2018                |
| <b>Balance:</b>  | Sartorius -<br>LE244S-OCE |

