

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 715 010

**Graduated cylinder, 100:1,0 ml, A, conf. cert.**




**Batch:** 25-18

**Serial:**

**Measuring Nr.** 21841

**Batch size:** 2400

Volume:	20,000 ml	60,000 ml	100,000 ml
<b>Tolerance:</b>	± 0,5000 ml	± 0,5000 ml	± 0,5000 ml
<b>Average value:</b>	19,9514 ml	59,9812 ml	100,0250 ml
<b>Deviation:</b>	0,1068 ml	0,0780 ml	0,0691 ml
<b>Precision (CV):</b>	0,5342 %	0,1300 %	0,0691 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

<b>Operator:</b>	4
<b>Date:</b>	27.06.2018
<b>Balance:</b>	Sartorius - LA 1200 S

