

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 010

Graduated cylinder, 100:1,0 ml, A, conf. cert.




Batch: 37-20

Serial:

Measuring Nr.: 40320

Batch size: 1200

| Volume: | 20,000 ml | 60,000 ml | 100,000 ml |
|------------------------|-------------|-------------|-------------|
| Tolerance: | ± 0,5000 ml | ± 0,5000 ml | ± 0,5000 ml |
| Average value: | 20,0136 ml | 60,0318 ml | 100,0446 ml |
| Deviation: | 0,0455 ml | 0,0529 ml | 0,0714 ml |
| Precision (CV): | 0,2274 % | 0,0881 % | 0,0714 % |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

| | |
|------------------|---------------------------|
| Operator: | 7 |
| Date: | 22.09.2020 |
| Balance: | Sartorius - LE244S-OCE |

