

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 010

Graduated cylinder, 100:1,0 ml, A, conf. cert.




Batch: 32-21

Serial:

Measuring Nr.: 48385

Batch size: 1196

Volume:	20,0000 ml	60,0000 ml	100,0000 ml
Tolerance:	± 0,5000 ml	± 0,5000 ml	± 0,5000 ml
Average value:	19,9921 ml	59,9782 ml	99,9998 ml
Deviation:	0,1109 ml	0,0992 ml	0,0892 ml
Precision (CV):	0,5543 %	0,1653 %	0,0892 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	17.08.2021
Balance:	Sartorius - LA 1200 S

