

Certificate of performance



Article: 3 715 010

Graduated cylinder, 100:1,0 ml, A, conf. cert.




Batch: 12-24

Serial:

Measuring Nr. 66461

Batch size: 1298

Volume:	20,000 ml	60,000 ml	100,000 ml
Tolerance:	± 0,5000 ml	± 0,5000 ml	± 0,5000 ml
Average value:	20,0193 ml	60,0963 ml	100,0126 ml
Deviation:	0,1060 ml	0,0850 ml	0,0916 ml
Precision (CV):	0,5300 %	0,1417 %	0,0916 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	26.03.2024
Balance:	Sartorius - LE244S-OCE

