

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 025 BL

Grad. cylinder, 250:2 ml, blue grad., A, conf. cert.




Batch: 18-20

Serial:

Measuring Nr. 37563

Batch size: 500

| Volume: | 50,0000 ml | 150,0000 ml | 250,0000 ml |
|------------------------|-------------|-------------|-------------|
| Tolerance: | ± 1,0000 ml | ± 1,0000 ml | ± 1,0000 ml |
| Average value: | 49,8476 ml | 150,0445 ml | 250,1977 ml |
| Deviation: | 0,1765 ml | 0,1331 ml | 0,1001 ml |
| Precision (CV): | 0,3530 % | 0,0887 % | 0,0401 % |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

| | |
|------------------|--------------------------|
| Operator: | 4 |
| Date: | 29.04.2020 |
| Balance: | Sartorius - LA 1200 S |

