

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 025 BL

Grad. cylinder, 250:2 ml, blue grad., A, conf. cert.




Batch: 13-23

Serial:

Measuring Nr. 62394

Batch size: 279

Volume:	50,0000 ml	150,0000 ml	250,0000 ml
Tolerance:	± 1,0000 ml	± 1,0000 ml	± 1,0000 ml
Average value:	49,9208 ml	150,1089 ml	250,3459 ml
Deviation:	0,1268 ml	0,1147 ml	0,0999 ml
Precision (CV):	0,2537 %	0,0765 %	0,0400 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator: 7

Date: 04.04.2023

Balance: Sartorius -
LE244S-OCE

