

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 715 025 USP

**Graduated cylinders, 250 ml, A, conf. cert.**




**Batch:** 47-17

**Serial:**

**Measuring Nr.** 16172

**Batch size:** 50

Volume:	50,0000 ml	150,0000 ml	250,0000 ml
<b>Tolerance:</b>	± 1,0000 ml	± 1,0000 ml	± 1,0000 ml
<b>Average value:</b>	50,0549 ml	150,1096 ml	250,0442 ml
<b>Deviation:</b>	0,0985 ml	0,0902 ml	0,0434 ml
<b>Precision (CV):</b>	0,1971 %	0,0601 %	0,0174 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

<b>Operator:</b>	4
<b>Date:</b>	04.12.2017
<b>Balance:</b>	Sartorius - LA 1200 S

