

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 025 USP

Graduated cylinders, 250:2 ml, A, USP, conf. cert.

Batch: 42-20




Serial:

Measuring Nr.: 41591

Batch size:

100

Volume:	50,0000 ml	150,0000 ml	250,0000 ml
Tolerance:	± 1,0000 ml	± 1,0000 ml	± 1,0000 ml
Average value:	50,2045 ml	150,1676 ml	249,8269 ml
Deviation:	0,1147 ml	0,1844 ml	0,1562 ml
Precision (CV):	0,2294 %	0,1229 %	0,0625 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	21.10.2020
Balance:	Sartorius - LA 1200 S

