

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 715 025 USP

**Graduated cylinders, 250:2 ml, A, USP, conf. cert.**

**Batch:** 13-23




**Serial:**

**Measuring Nr.** 62408

**Batch size:**

280

Volume:	50,000 ml	150,000 ml	250,000 ml
<b>Tolerance:</b>	± 1,000 ml	± 1,000 ml	± 1,000 ml
<b>Average value:</b>	50,1756 ml	150,2177 ml	250,2675 ml
<b>Deviation:</b>	0,1565 ml	0,1271 ml	0,1332 ml
<b>Precision (CV):</b>	0,3131 %	0,0847 %	0,0533 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

**Operator:**

7

**Date:**

04.04.2023

**Balance:**

Sartorius -  
LE244S-OCE

