

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 025

Graduated cylinders, 250 ml, A, conf. cert.




Batch: 46-17

Serial: -

Measuring Nr.: 15979

Batch size: 800

Volume:	50,0000 ml	150,0000 ml	250,0000 ml
Tolerance:	± 1,0000 ml	± 1,0000 ml	± 1,0000 ml
Average value:	50,0397 ml	150,0392 ml	250,0331 ml
Deviation:	0,3038 ml	0,1772 ml	0,1972 ml
Precision (CV):	0,6077 %	0,1182 %	0,0789 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	28.11.2017
Balance:	Sartorius - LA 1200 S

