

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 025

Graduated cylinders, 250:2 ml, A, conf. cert.




Batch: 50-18

Serial:

Measuring Nr.: 26476

Batch size: 890

| Volume: | 50,000 ml | 150,000 ml | 250,000 ml |
|------------------------|------------|-------------|-------------|
| Tolerance: | ± 1,000 ml | ± 1,000 ml | ± 1,000 ml |
| Average value: | 49,8874 ml | 150,0024 ml | 250,0532 ml |
| Deviation: | 0,2186 ml | 0,1360 ml | 0,1573 ml |
| Precision (CV): | 0,4372 % | 0,0907 % | 0,0629 % |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

| | |
|------------------|--------------------------|
| Operator: | 4 |
| Date: | 13.12.2018 |
| Balance: | Sartorius - LA 1200 S |

