

Certificate of performance



Article: 3 715 025

Graduated cylinder, 250:2 ml, A, conf. cert.




Batch: 43-23

Serial:

Measuring Nr. 64915

Batch size: 280

Volume:	50,000 ml	150,000 ml	250,000 ml
Tolerance:	± 1,000 ml	± 1,000 ml	± 1,000 ml
Average value:	50,1447 ml	150,1423 ml	250,1597 ml
Deviation:	0,2256 ml	0,2228 ml	0,2266 ml
Precision (CV):	0,4512 %	0,1486 %	0,0906 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	07.11.2023
Balance:	Sartorius - LA 1200 S

