

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 050 BL USP

Grad. cylinder, 500:5 ml, blue grad., USP, A, conf. cert.




Batch: 06-19

Serial: -

Measuring Nr.: 27567

Batch size: 100

Volume:	100,000 ml	300,000 ml	500,000 ml
Tolerance:	± 2,0000 ml	± 2,0000 ml	± 2,0000 ml
Average value:	99,6092 ml	300,0009 ml	499,9948 ml
Deviation:	0,2482 ml	0,1750 ml	0,2524 ml
Precision (CV):	0,2482 %	0,0583 %	0,0505 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

4

Date:

07.02.2019

Balance:

Sartorius - LA
1200 S

DE-M