

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 050 BL USP

Grad. cylinder, 500:5 ml, blue grad., USP, A, conf. cert.




Batch: 30-21

Serial:

Measuring Nr.: 48074

Batch size: 150

Volume:	100,000 ml	300,000 ml	500,000 ml
Tolerance:	± 2,0000 ml	± 2,0000 ml	± 2,0000 ml
Average value:	99,4786 ml	299,7927 ml	500,0343 ml
Deviation:	0,3767 ml	0,4131 ml	0,2042 ml
Precision (CV):	0,3767 %	0,1377 %	0,0408 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	27.07.2021
Balance:	Sartorius - LE244S-OCE

