

Certificate of performance



Article: 3 715 050 BL

Grad. cylinders, 500:5 ml, blue grad., A, conf. cert.




Batch: 20-18

Serial:

Measuring Nr. 20609

Batch size: 280

Volume:	100,000 ml	300,000 ml	500,000 ml
Tolerance:	± 2,5000 ml	± 2,5000 ml	± 2,5000 ml
Average value:	100,0936 ml	300,3486 ml	500,6169 ml
Deviation:	0,4346 ml	0,3401 ml	0,3194 ml
Precision (CV):	0,4346 %	0,1134 %	0,0639 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	44
Date:	15.05.2018
Balance:	Sartorius - LE244S-OCE

