

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 050 BL

Grad. cylinder, 500:5 ml, blue grad., A, conf. cert.




Batch: 20-23

Serial:

Measuring Nr. 63138

Batch size: 150

Volume:	100,000 ml	300,000 ml	500,000 ml
Tolerance:	± 2,5000 ml	± 2,5000 ml	± 2,5000 ml
Average value:	100,4310 ml	299,9246 ml	499,6090 ml
Deviation:	0,5505 ml	0,2912 ml	0,2375 ml
Precision (CV):	0,5505 %	0,0971 %	0,0475 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	23.05.2023
Balance:	Sartorius - LE244S-OCE

