

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 050 USP

Grad.cylinder, 500:5 ml, USP, A, conf. cert.




Batch: 42-19

Serial:

Measuring Nr. 33743

Batch size: 80

Volume:	100,000 ml	300,000 ml	500,000 ml
Tolerance:	± 2,0000 ml	± 2,0000 ml	± 2,0000 ml
Average value:	100,0917 ml	300,1717 ml	500,0896 ml
Deviation:	0,0637 ml	0,5009 ml	0,4387 ml
Precision (CV):	0,0637 %	0,1670 %	0,0877 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	17.10.2019
Balance:	Sartorius - LA 1200 S

DE-M