

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 050 USP

Grad. cylinder, 500:5 ml, USP, A, conf. cert.




Batch: 11-20

Serial: -

Measuring Nr.: 36658

Batch size: 100

Volume:	100,000 ml	300,000 ml	500,000 ml
Tolerance:	± 2,0000 ml	± 2,0000 ml	± 2,0000 ml
Average value:	99,9371 ml	300,2035 ml	500,2221 ml
Deviation:	0,4730 ml	0,3695 ml	0,4653 ml
Precision (CV):	0,4730 %	0,1232 %	0,0931 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

4

Date:

13.03.2020

Balance:

Sartorius - LA
1200 S

DE-M