

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 715 050 USP

**Grad. cylinder, 500:5 ml, USP, A, conf. cert.**




**Batch:** 30-21

**Serial:**

**Measuring Nr.:** 48075

**Batch size:** 280

Volume:	100,000 ml	300,000 ml	500,000 ml
<b>Tolerance:</b>	± 2,0000 ml	± 2,0000 ml	± 2,0000 ml
<b>Average value:</b>	100,0975 ml	300,4568 ml	500,3810 ml
<b>Deviation:</b>	0,2366 ml	0,2591 ml	0,3226 ml
<b>Precision (CV):</b>	0,2366 %	0,0864 %	0,0645 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

<b>Operator:</b>	7
<b>Date:</b>	28.07.2021
<b>Balance:</b>	Sartorius - LA 1200 S

