

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 050 USP

Grad. cylinder, 500:5 ml, USP, A, conf. cert.




Batch: 44-22

Serial:

Measuring Nr. 60185

Batch size: 280

Volume:	100,000 ml	300,000 ml	500,000 ml
Tolerance:	± 2,0000 ml	± 2,0000 ml	± 2,0000 ml
Average value:	99,8895 ml	300,1942 ml	500,0756 ml
Deviation:	0,3356 ml	0,2957 ml	0,1534 ml
Precision (CV):	0,3356 %	0,0986 %	0,0307 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

4

Date:

08.11.2022

Balance:

Sartorius - LA
1200 S

DE-M