

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 715 050 USP

**Grad. cylinder, 500:5 ml, USP, A, conf. cert.**




**Batch:** 20-23

**Serial:**

**Measuring Nr.:** 63139

**Batch size:** 100

| Volume:         | 100,000 ml  | 300,000 ml  | 500,000 ml  |
|-----------------|-------------|-------------|-------------|
| Tolerance:      | ± 2,0000 ml | ± 2,0000 ml | ± 2,0000 ml |
| Average value:  | 100,4114 ml | 299,9080 ml | 499,8168 ml |
| Deviation:      | 0,3531 ml   | 0,3031 ml   | 0,1949 ml   |
| Precision (CV): | 0,3531 %    | 0,1010 %    | 0,0390 %    |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

|           |                           |
|-----------|---------------------------|
| Operator: | 7                         |
| Date:     | 24.05.2023                |
| Balance:  | Sartorius -<br>LE244S-OCE |

