

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 050

Grad. cylinder 500 ml class A




Batch: 33-17

Serial:

Measuring Nr.: 13895

Batch size: 500

Volume:	100,000 ml	300,000 ml	500,000 ml
Tolerance:	± 2,5000 ml	± 2,5000 ml	± 2,5000 ml
Average value:	100,2540 ml	300,1006 ml	499,9667 ml
Deviation:	0,5023 ml	0,2431 ml	0,3807 ml
Precision (CV):	0,5023 %	0,0810 %	0,0761 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	06.09.2017
Balance:	Sartorius - LA 1200 S

